



Glo-Tip[®]

SPRAY CATHETER



Deliver a fine staining-solution mist.

When spraying the gastrointestinal mucosa to identify suspicious lesions, a fine mist is essential to optimally stain the tissue. The Glo-Tip Spray Catheter delivers an exceptionally fine mist, allowing you to perform effective procedures. The 7 FR Teflon[®] sheath provides low coefficient of friction for even flow of the staining solution and for ease of insertion.

www.cookmedical.com



Glo-Tip®

SPRAY CATHETER

This device is used in conjunction with an endoscope to spray the gastrointestinal mucosa.

Supplied sterile and is disposable - intended for single use only.

- 7 FR Teflon® sheath provides low coefficient to friction for ease of insertion, allows for even flow of staining solution
- Radiopaque tip enhances fluoroscopic visualization, if needed
- Smooth metal spray tip is atraumatic to mucosa
- Disposability eliminates cleaning and reprocessing costs, saves time

Order Number	GPN	French Size	Length	Tip
GT-7-SPRAY	G24892	7 Fr	240 cm	Spray

Color Code Yellow - This device will not accept a wire guide.



COOK ENDOSCOPY

4900 Bethania Station Road, Winston-Salem, NC 27105 U.S.A.
Phone: 336 744-0157, Toll Free: (USA) 800 457-4500, Fax: 336 744-5231

Glo-Tip is a registered trademark of Cook Endoscopy.
Cook is a registered trademark of Cook Incorporated.

© 2010 Cook Medical

L-17834/0310

www.cookmedical.com

AORTIC
INTERVENTION

CRITICAL
CARE

ENDOSCOPY

INTERVENTIONAL
RADIOLOGY

LEAD
MANAGEMENT

PERIPHERAL
INTERVENTION

SURGERY

UROLOGY

WOMEN'S
HEALTH